# PK macros translation rules (with examples)

Maciej Swat (EMBL-EBI) 16 April 2015

### **Assumptions**

All translation rules are based on the following assumptions:

- 1. the PK macros to be translated are correctly encoded this means that no validation rules are formulated and enforced (for now) and that the validation should be defined as a separate task
- 2. no defaults are used, i.e.
  - a. macros use always the 'cmt' argument, where applicable
  - b. (a) means that following macro is valid: **iv()** which corresponds to **iv(cmt=1)** it is when only one compartment is defined and therefore the 'cmt' argument is not required
  - c. 'adm/type' is always required even if there is a single route of administration
  - d. 'amount' argument in the <peripheral> macro is required

## Proposed strategy

- 1. Start with real-life examples, so called NONMEM PREDPP library models: ADVAN1-4 & 10-12 (see example section below). This is a set of well-known and frequently used models in pharmacometrics (described in detail in the PK macro report [2]), then go to more complex examples.
- 2. Models mentioned in (1) must be correctly formulated as full macros, i.e. no defaults.
- 3. When testing more general macros, making use of defaults, one should
  - a. first, complete the given set of macros with missing arguments, e.g. instead of 'iv()' write 'iv(cmt=1)' etc.
  - b. process then the complete macros using the translation rules.

### Translation rules

This is a set of rules collected while doing the R converter for the PK macros [3]. See [2] for the full set of macros and arguments. The original document [1] from Lixoft describes the elements and their meaning in detail.

- 1. Start with (1) 'compartment' and (2) 'peripheral' macros in that order they define the core structure followed by remaining ones
- 2. Identify associated
  - a. compartment numbers, e.g. 'cmt=1' store in cmtNumber array
  - b. amounts names, e.g. 'amount=Ac' store in cmtAmount array
  - c. volume names, e.g. 'volume=V' store in cmtVolume array
  - d. concentration names, e.g. 'concentration=C' store in cmtConcentration array
- 3. For each 'compartment' macro, e.g.

```
compartment(cmt=1, amount=Ac, volume=V, concentration=C)
create an empty ODE 'dcmtAmount[cmt]/dt = ', i.e. 'dAc/dt = '
```

4. For each 'peripheral' macro, e.g.

- a. create an empty 'j' ODE 'dcmtAmount[j]/dt = ', i.e. 'dAp/dt = '
- b. process 'kij' or 'k\_i\_j'arguments label 'i' is for one of the central 'compartment's defined before, label 'j' is for the current peripheral compartment
  - i. Add '- kij \* cmtAmount[j]' to the 'i' ODE
  - ii. Add '+ kij \* cmtAmount[i]' to the 'current' ODE
- c. process 'kji' or 'k\_j\_i'arguments
  - i. Add '- kji \* cmtAmount[i]' to the 'j' ODE
  - ii. Add '+ kji \* cmtAmount[j]' to the 'current' ODE
- 5. 'absorption/oral' macro three options: {either 'Tk0', 'ka' or 'ka, Ktr,Mtt'}

- a. each such macro means a new 'depot' compartment
- b. create new ODE string 'dcmtAmount[new name]/dt = ', i.e. 'dAd/dt = '
- c. identify the target compartment, 'i'
- d. update compartment/amount 'arrays'
- e. case1: zero order absorption, example 9

- i. add '- Tk0' to the new ODE string
- ii. add '+ Tk0' to the 'i' target compartment ODE
- f. case2: first order absorption must have 'ka' argument

- i. add '+ ka\*cmtAmount[new 'j']' to the target compartment 'i'
- ii. add '- ka\*cmtAmount[new 'j']' to the current depot compartment ODE
- g. case3: (example 11) models with *transit* compartments absorption have additional 'Ktr' and 'Mtt'

- i. new absorption compartment, e.g. 'Aa'
- ii. add '+ ka\*Aa' to the target compartment 'i'
- iii. add the following new ODE 'dAa/dt = exp[log(F\*Dose)) + log(Ktr) + n\*log(Ktr\*(t-t\_Dose)) Ktr\*(t-t\_Dose) log(n!)] ka\*Aa'
- iv. for simplicity (ii) assumes there is only one administration defined with 'Ktr/Mtt', otherwise need more As's compartments and related ODs, i.e. Aa1 with 'dAa1/dt=...'
- v. target 'Dose'
- h. As the last step one needs to provide the information about the input, required for the link with a dataset
  - i. the input is 'oral'
  - ii. administration number is given by 'adm/type'
  - iii. target cmtAmount[i] (except case 3)
  - iv. i-iii: New Input[inputNumber] ORAL administration, adm=a, target=cmtAmount[i]
- 6. 'iv' macro

- a. in this case there are NO additions to ODEs
- b. only the input information needs to be provided, similar to the 'oral' case
  - i. the input is 'iv'
  - ii. administration number is given by 'adm/type'
  - iii. target cmtAmount[i]
  - iv. i-iii: Input[inputNumber] IV administration, adm=a, target=cmtAmount[i]

7. 'transfer' macro

- a. note here an assignment for kt, 'kl' to be used
- b. extract 'to' target compartment number, 'from' source compartment number
- c. assign to 'to' compartment ' + kl\*cmtAmount[i]'
- d. assign to 'from' compartment ' kl\*cmtAmount[j]'
- 8. 'elimination' macro three options {either 'k', 'Km&Vm' or 'CL&V'}
  - a. case 1: 'k' linear elimination

add to the 'cmt' compartment the ' - k\*cmtAmount[i]'

b. case 2: linear elimination with CL

add to the 'cmt' compartment the ' - CL/V\*cmtAmount[i]'

c. case3: 'Km & Vm' saturable elimination

add to the 'cmt' compartment the ' - Vm\*cmtAmount[i]/(Km + cmtAmount[i])'

9. 'effect' macro - see example 8

- a. update compartment/volume/amount 'arrays' although the last two will have 'NaN's for this macro
- b. create new algebraic equation cmtConcentration[i] = cmtAmount[i]/cmtVolume[i],i.e. C= Ac/V
- c. create a new ODE 'dCe/dt = ke0\*(cmtConcentration[i] Ce)'
- 10. '**depot**' macro is a bit special it occurs only in connection with explicitly defined ODE's in the MLXTRAN literature, e.g.

**PK:** depot(adm=a, target=Ac)

**EQUATION:** ddt Ac = -k\*Ac

and means bolus IV administration

or

PK: depot(adm=a, target=Ac, ka)

**EQUATION**: ddt Ac = -k\*Ac

and means ORAL administration

a. case 1: without 'ka' argument

no additions to the ODEs - only to 'Input[inputNumber]' for a new 'inputNumber'

- i. New Input[inputNumber] IV administration, adm=a, target=cmtAmount[i]
- b. case 2: with 'ka' argument

depot(adm=a, target=Ac, ka)

corresponds to these 2 macros

- i. i.e. creates new *depot* compartment and according ODE 'd*cmtAmount*[new depot name]/dt = ', i.e. 'dAd/dt = ka\*Ad'
- ii. adds ' + ka\*Ad' to the target compartment, Ac.
- iii. NEW Input[inputNumber] ORAL administration, adm=a, target=cmtAmount[i]

### **Examples**

The following 'ODEs' and 'Input's have been generated by the set of R scripts - see [3] - NOTE, the scripts underwent a number of changes - to be updated soon.

### Example 1: ADVAN1, TRANS1

Macros: compartment(cmt=1, amount=Ac, volume=V)

iv(adm=1, cmt=1)
elimination(cmt=1, k)

**ODEs**: dAc/dt = - k\*Ac

Input: Input[1]: IV administration, adm=1, target=Ac

### Example 2: ADVAN2, TRANS1

Macros: compartment(cmt=1, amount=Ac, volume=V)

oral(adm=1, cmt=1, ka)
elimination(cmt=1, k)

**ODEs:** dAc/dt = + ka\*Ad2 - k\*Ac

dAd2/dt = - ka\*Ad2

Input: Input[1]: ORAL administration, adm=1, target=Ad2

### Example 3: ADVAN3, TRANS1

Macros: compartment(cmt=1, amount=Ac, volume=V)

peripheral(k12, k21, amount=Ap)

iv(adm=1, cmt=1)
elimination(cmt=1, k)

**ODEs**: dAc/dt = - k12\*Ac + k21\*Ap - k\*Ac

dAp/dt = k12\*Ac - k21\*Ap

Input:
Input[1]: IV administration, adm=1, target=Ac

#### Example 4: ADVAN4, TRANS1

Macros: compartment(cmt=1, amount=Ac, volume=V)

peripheral(k12, k21, amount=Ap)

oral(adm=1, cmt=1, ka)
elimination(cmt=1, k)

**ODEs:** dAc/dt = - k12\*Ac + k21\*Ap + ka\*Ad3 - k\*Ac

dAp/dt = k12\*Ac - k21\*Ap

dAd3/dt = - ka\*Ad3

Input: Input[1]: ORAL administration, adm=1, target=Ad3

## Example 5: ADVAN10, TRANS1

Macros: compartment(cmt=1, amount=Ac, volume=V)

```
iv(adm=1, cmt=1)
elimination(cmt=1, Km, Vm)
```

ODEs: dAc/dt = - Vm\*Ac/(Km + Ac)

Input:
Input[1]: IV administration, adm=1, target=Ac

### Example 6: ADVAN11, TRANS1

Macros: compartment(cmt=1, amount=Ac, volume=V)

peripheral(k12, k21, amount=Ap1)
peripheral(k13, k31, amount=Ap2)

iv(adm=1, cmt=1)
elimination(cmt=1, k)

ODEs: dAc/dt= - k12\*Ac + k21\*Ap1 - k13\*Ac + k31\*Ap2 - k\*Ac

dAp1/dt= k12\*Ac - k21\*Ap1 dAp2/dt= k13\*Ac - k31\*Ap2

Input:
Input[1]: IV administration, adm=1, target=Ac

### Example 7: ADVAN12, TRANS1

Macros: compartment(cmt=1, amount=Ac, volume=V)

peripheral(k12, k21, amount=Ap1)
peripheral(k13, k31, amount=Ap2)

oral(adm=1, cmt=1, ka)
elimination(cmt=1, k)

**ODEs:** dAc/dt = -k12\*Ac + k21\*Ap1 - k13\*Ac + k31\*Ap2 + ka\*Ad4 - k\*Ac

dAp1/dt= k12\*Ac - k21\*Ap1 dAp2/dt= k13\*Ac - k31\*Ap2

dAd4/dt = - ka\*Ad4

Input: Input[1]: ORAL administration, adm=1, target=Ad4

#### Example 8: Model with effect compartment

Macros: compartment(cmt=1, amount=Ac, volume=V, concentration=C)

iv(adm=1, cmt=1)
elimination(cmt=1, k)

effect(cmt=1, ke0, concentration=Ce)

**ODEs**: dAc/dt = - k\*Ac

dCe/dt = ke0\*(C - Ce)

C=Ac/V

Input:
Input[1]: IV administration, adm=1, target=Ac

## Example 9: Model with oral and Tk0

Macros: compartment(cmt=1, amount=Ac, concentration=Cc, volume=V)

oral(adm=1, cmt=1, Tk0)
elimination(cmt=1, k)

ODEs: dAc/dt = + Tk01 - k\*Ac

dAd2/dt = - Tk01

Input: Input[1]: ORAL administration, adm=1, target=Ad2

### Example 10: Sequential zero order/first order absorption processes - only one 'adm'

Macros: compartment(cmt=1, amount=Ac, concentration=Cc, volume=V)

oral(adm=1, cmt=1, Tk0, p=F0)

oral(adm=2, cmt=1, ka , Tlag=Tk0 , p=1-F0)

elimination(cmt=1, k)

**ODEs:** dAc/dt = + Tk01 + ka\*Ad3 - k\*Ac

dAd2/dt = - Tk01dAd3/dt = - ka\*Ad3

Input: Input[1]: ORAL administration, adm=1, target=Ad1; p=F0

Input[2]: ORAL administration, adm=2, target=Ad3; Tlag=Tk0; p=1-F0

### Example 11: Model with transit compartments. 'example\_1comp\_kaKtrMtt\_k.txt'

Macros: compartment(cmt=1,amount=Ac,volume=V,concentration=C)

oral(adm=1, cmt=1, Mtt, Ktr, ka)

elimination(cmt=1, k)

**ODEs**: dAc/dt = + ka\*Aa - k\*Ac

dAa/dt = exp[log(F\*Dose)) + log(Ktr) + n\*log(Ktr\*(t-t Dose)) - Ktr\*(t-t Dose)

 $-\log(n!)$ ] -ka\*Aa

Input:
Input[1]: ORAL administration, adm=1, target=Dose

## References

- [1] MLXTRAN\_forMonolix\_May2014.pdf copy in the folder 'PK macros'
- [2] PKmacros\_in\_PharmML0.6\_17Feb2015.pdf copy in the folder 'PK macros'
- [3] PKmacros2ODEs.zip R converter copy in the folder 'PK macros'
  - run the 'testRun.R' script, the macro sets for ADVAN models are defined in line 10.